

Rinse

Acidic Neutralizing Rinse

ID 4103

Description

For use in automated cleaning of laboratory glassware in medical and research facilities. Neutralizes residual alkalinity following wash function. Can also be used as acidic detergent for washing stainless animal cages. Helps to remove lime accumulations in automated washers. Do not use on anodized aluminum or light metal alloys.

Directions

Use at a rate of 1-2 mls per gallon of rinse water/tank capacity.

Technical Data

Form.....	Liquid	Rinsability	Excellent
Color	Clear	Solubility	Rapid
Scent.....	Mild Acidic	pH (1% solution)	~2.8



WARNING

HAZARDS: Causes serious eye irritation.

PRECAUTIONS: Wash hands thoroughly after handling. Wear protective gloves and eye protection.

RESPONSE/FIRST AID: Have product container or label with you when calling poison control center or physician. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, get medical advice/attention. **IF ON SKIN:** Wash with plenty of soap and water. Take off contaminated clothing and wash before re-use. If skin irritation occurs, get medical advice/attention. **IF SWALLOWED:** Rinse mouth with water. Call a poison control center or physician for advice.

STORAGE & DISPOSAL: No GHS phrases apply. Dispose of contents/container in accordance with all local, regional, national and international regulations.

