

Treat

Enzymatic Pretreat

ID 4688

Description

Treat is the ideal surgical instrument enzymatic pre-cleaner, prevents bioburden build up. Perfect for Operating Rooms, Central Supply, Endoscopy Suites and Outpatient Surgery. The Treat enhanced foaming system produces super clinging foam for outstanding soil penetration.

Directions

Spray Treat onto instruments immediately following operation. Let instruments soak before & during transport to cleaning area. Apply additional Treat as needed. Instruments must be followed up with normal cleaning and sterilization procedures as Treat is not a sterilizer.

Technical Data

Form.....	Liquid	Scent.....	Mild Surfactant
Color	Clear, blue, viscous	pH (1% solution)	8-9

HAZARDS & PRECAUTIONS: No GHS phrases apply.

RESPONSE/FIRST AID: PROLONG SKIN CONTACT: Wash skin with clear water and mild soap for 15 minutes. Contact physician if irritation occurs. **INDIGESTION:** Wash out mouth with water provided person is conscious. Seek medical attention if feeling unwell. **EYE CONTACT:** Flush with water for 15 minutes. If irritation persists, get medical attention.

STORAGE & DISPOSAL: No GHS phrases apply. Dispose of contents/container in accordance with all local, regional, national and international regulations.



Manufactured for:

Spire Integrated Solutions

A Chalmers Group Company

8719 S 135th St, Ste #300, Omaha, NE 68135 | (402) 344-4200